## Anaphylaxis Management: School Supervisors' Observation Checklist



An observation record must be made and retained at the school for <u>each</u> staff member demonstrating the correct use of the adrenaline autoinjector (trainer) device. Certification that the device is used correctly can only be provided by the appropriately trained School Anaphylaxis Supervisor.

Name of School Anaphylaxis Supervisor:	Signature:
Name of staff member being assessed:	Signature:
Assessment Result:	Competent or Not competent (select as appropriate)
Assessment date:	

## Verifying the correct use of Adrenaline Autoinjector (trainer) Devices

Stage	Actions	√orx
Preparation	Verification resources, documentation and adrenaline autoinjector (trainer) devices and equipment are on hand and a suitable space for verification is identified.	
Preparation	Confirmation of the availability of a mock casualty (adult) for the staff member to demonstrate use of the adrenaline autoinjector devices on.  Testing of the device on oneself or the verifier is not appropriate.	
Demonstration	Successful completion of the ASCIA Anaphylaxis e-training for Victorian Schools within the previous 30 days is confirmed by sighting the staff member's printed ASCIA e-training certificate.	
Demonstration	Confirmation the staff member has had an opportunity to practise use of the adrenaline autoinjector (trainer) device/s prior to the verification stage.	
Demonstration	To conduct a fair appraisal of performance, the verifier should first explain what the candidate is required to do and what they will be observed doing prior to the demonstration, including a scenario for the mock casualty. This ensures the candidate is ready to be verified and clearly understands what constitutes successful performance or not.	

## **Practical Demonstration**

Stage	The staff member:	Attempts ✓ or ×		
		1	2	3
Prior to use:  Identifying the components of the EpiPen®	Correctly identified components of the adrenaline autoinjector (although some of these are not available on the trainer device, they should be raised and tested): School Anaphylaxis Supervisors to ask the below questions.  • Where is the needle located?  • What is a safety mechanism of the EpiPen®?  • What triggers the EpiPen® to administer the medication?  • What does the label of the EpiPen® show?			

Stage	The staff member:	Atten	Attempts ✓ or ×	
		1	2	3
Prior to use:  Demonstrated knowledge of the appropriate checks of the EpiPen®	Demonstrated knowledge of the appropriate checks of the adrenaline autoinjector device (although these are not available on the trainer device, they should be raised and tested): School Anaphylaxis Supervisors to ask the below questions.  • Prior to administering the EpiPen® what should you check?  • What do you check the viewing window for?  • What do you check the label for?			
Correct positioning when applying anaphylaxis first aid	Positioned themselves and the (mock) casualty correctly in accordance with ASCIA guidelines ensuring the:  casualty is lying flat unless breathing is difficult or placed in a recovery position if unconscious or vomiting  casualty is securely positioned to prevent movement when administering the adrenaline autoinjector device  person administering the adrenaline autoinjector device is facing the casualty.			
Demonstration:  Correct administration of the EpiPen®	Administered the adrenaline autoinjector device correctly (this example is for an EpiPen® device):  • formed a fist to hold the EpiPen® device correctly  • pulled off blue safety release  • applied the orange end at right angle to the outer mid-thigh (with or without clothing), ensuring pockets and seams were not in the way  • activated the EpiPen® by pushing down hard until a click is heard  • hold the EpiPen® in position for 3 seconds after activation  • removed EpiPen®.			
Demonstration:	Demonstrated correct use in a realistic time period for treatment in an emergency situation.			
Post use:  Handling used  EpiPen®	Demonstrated knowledge of correct procedures post use of the adrenaline autoinjector devices: School Anaphylaxis Supervisors to ask the below questions.  • What information should you record at the time of administering the EpiPen®?  • What do you do with the used EpiPen® once it has been			

## **Test Outcome**

	Certifying the correct use of the adrenaline autoinjector (training) device	√or×
	Where checking and demonstration is successful the verifier will:	
ı	sign and date the staff member's ASCIA e-training certificate	
	provide a copy to the staff member	
	store the staff member's ASCIA certificate and this observation record in a central office location to	
ı	ensure confidentiality is maintained	
ı	update school staff records for anaphylaxis training.	
	If the adrenaline autoinjector (trainer) device has <b>NOT</b> been checked or administered correctly through	
ı	successfully completing all the steps above, the verifier cannot deem the staff member competent. The	
1	staff member should be referred back to the ASCIA Anaphylaxis e-training for further training and re-	
ı	present for verification:	
	this action should be recorded in staff records	
ı	the verifier must not provide training to correct practice.	